

LOVIBOND® WATER ANALYSIS

TINTOMETER® GROUP



Lovibond® ShuddhCheck™

Your gateway to better Water Testing

INSTRUCTION GUIDE- pH COLOUR CARD

WORKING RANGE: 2 – 13

MEASURING RANGE: 2, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13



Test Procedure

1. Take 5 ml water sample in a glass vial
2. Add 1 drop of pH Indicator and mix well
3. Wait for 5 minutes
4. Notice the colour change and compare it with the colour card provided with the kit
5. The matching colour corresponds to the pH of the water sample

This will give the exact pH value directly of the water sample.

No. of Tests: 100

Order Code: IN415020MK

Reorder Code: IN415021MR

Tintometer India Private Limited

Door No. 7-2-C-14, 2nd, 3rd & 4th Floor

Sanathnagar Industrial Estate | Hyderabad-500018

Telangana | India



Notes

- * Cleaning: Thoroughly rinse out sample container after use.
- ** Test reagents are for chemical testing only. Not meant for human consumption; Keep out of reach of children. For industrial use only.



For more information please reach out to us:

☎ 1800 599 3891/92 | ✉ indiaoffice@lovibond.in | 🌐 www.lovibond.in

🌐 <https://in.linkedin.com/company/tintometer-india-private-limited>